

Your SureMed MultiMed 31-Day Starter Kit



Your SureMed 31-day starter kit includes: 26 cards, 26 blisters, a filling foam and sealing roller.

Tips


- This blister packaging is not child resistant, please store beyond the reach of children.
- Keep unused cards in a clean, dry and air conditioned location.
- If there is a slight tug while peeling off the pressure sensitive liner try peeling from a different corner.

Our Medications



We'll be filling 2 cards, 1 card with five morning medications and 1 card with four evening medications.

1. Remove Liner



Start with the morning card, remove the pressure sensitive liner revealing the sealing surface.

2. Filling Foam



Place the card on the filling foam template (dry glue side up).

3. Insert Blisters



Insert the blister sheet into the card.

4. Filling Morning - Pill #1



You take 5 different medications in the morning. Place this first medication in each of the 31 blisters.

Filling Morning - Pill #2



Place your second medication into each of the 31 blisters.

Filling Pill #3



You also take pill #3 every morning. Place one pill in each of the morning blisters, sunday thru saturday.

Filling Pill #4



Next, you take pill #4 at noon each day. Place one pill in each of the noon blisters, sunday thru saturday.

Filling Pill #5



You also take pill #5 each day at noon. Place one pill in each of the noon blisters, sunday thru saturday.

Filling Pill #6



Pill #6 you take each day before you go to bed. Place one pill in each of the bedtime blisters, sunday thru saturday.

Filling Pill #7



You also take pill #7 at bedtime each evening as well. Place one pill in each of the bedtime blisters, sunday thru saturday.

Fold



Once all of your medications are loaded into the card, fold the foil flap over to seal the card.

Seal



Using the roller, apply even pressure to seal.

SureMed 7-Day Finished



You're Done!
Your fresh clean disposable SureMed 7-Day card is filled and ready to use.

Qube Mini 7-Day Finished



You're Done!
Your fresh clean disposable Qube Mini 7-Day card is filled and ready to use.