

## Your SureMed MultiMed 31-Day Starter Kit



Your SureMed 31-day starter kit includes: 26 cards, 26 blisters, a filling foam and sealing roller.

## Tips


- This blister packaging is not child resistant, please store beyond the reach of children.
- Keep unused cards in a clean, dry and air conditioned location.
- If there is a slight tug while peeling off the pressure sensitive liner try peeling from a different corner.

## Our Medications



We'll be filling 2 cards, 1 card with five morning medications and 1 card with four evening medications.

## 1. Remove Liner



Start with the morning card, remove the pressure sensitive liner revealing the sealing surface.

## 2. Filling Foam



Place the card on the filling foam template (dry glue side up).

## 3. Insert Blisters



Insert the blister sheet into the card.

## 4. Filling Morning - Pill #1



You take 5 different medications in the morning. Place this first medication in each of the 31 blisters.

## Filling Morning - Pill #2



Place your second medication into each of the 31 blisters.



### Filling Pill #3



You also take pill #3 every morning. Place one pill in each of the morning blisters, sunday thru saturday.

### Filling Pill #4



Next, you take pill #4 at noon each day. Place one pill in each of the noon blisters, sunday thru saturday.

### Filling Pill #5



You also take pill #5 each day at noon. Place one pill in each of the noon blisters, sunday thru saturday.

### Filling Pill #6



Pill #6 you take each day before you go to bed. Place one pill in each of the bedtime blisters, sunday thru saturday.

### Filling Pill #7




You also take pill #7 at bedtime each evening as well. Place one pill in each of the bedtime blisters, sunday thru saturday.

### Fold




Once all of your medications are loaded into the card, fold the foil flap over to seal the card.

### Seal



Using the roller, apply even pressure to seal.

### SureMed 7-Day Finished



**You're Done!**  
Your fresh clean disposable SureMed 7-Day card is filled and ready to use.

### Qube Mini 7-Day Finished



**You're Done!**  
Your fresh clean disposable Qube Mini 7-Day card is filled and ready to use.