

**Your SureMed 7-Day or Qube Mini 7-Day Multi-Med Starter Kit**



SureMed 7-Day, Medicine-On-Time 7-Day & Qube Mini 7-Day kits includes: 26 blisters & cards, filling foam and sealing roller.

**Tips**

- Using the reverse side day and time of day as a guide fill the appropriate blister with the correct medication.
- This blister packaging is not child resistant, please store beyond the reach of children.
- Keep unused cards in a clean, dry and air conditioned location.
- If there is a slight tug while peeling off the pressure sensitive liner try peeling from a different corner.

**1. Remove Liner**



Remove the pressure sensitive liner revealing the sealing surface.  
**NOTE: Qube Card has a pressure liner on the middle panel only.**

**2. Filling Foam**



Place the card on the filling foam template (dry glue side up).

**3. Insert Blisters**



Insert the blister sheet into the card starting with the top row. Lower the rest of the sheet into the remaining holes.

**4. Filling Pill #1**



Say you take 7 medications a day and one pill twice a day. Pill #1 you take in the morning. Place one pill in each of the morning blisters, Sunday thru Saturday.

**Filling Pill #1**



You also take the same (pill #1) each evening as well. Place one pill in each of the evening blisters, Sunday thru Saturday.

**Filling Pill #2**



Pill #2 you take each morning. Place one pill in each of the morning blisters, Sunday thru Saturday.

### Filling Pill #3



You also take pill #3 every morning. Place one pill in each of the morning blisters, sunday thru saturday.

### Filling Pill #4



Next, you take pill #4 at noon each day. Place one pill in each of the noon blisters, sunday thru saturday.

### Filling Pill #5



You also take pill #5 each day at noon. Place one pill in each of the noon blisters, sunday thru saturday.

### Filling Pill #6



Pill #6 you take each day before you go to bed. Place one pill in each of the bedtime blisters, sunday thru saturday.

### Filling Pill #7



You also take pill #7 at bedtime each evening as well. Place one pill in each of the bedtime blisters, sunday thru saturday.

### Fold



Once all of your medications are loaded into the card, fold the foil flap over to seal the card.

### Seal



Using the roller, apply even pressure to seal.

### Our 3 Finished 7-Day 4-Pass Cards



Medicine-On-Time Card



SureMed Card



Qube Mini Card

### You're Done!

Your fresh clean disposable Medicine-On-Time 7-Day 4-Pass Card, the SureMed 7-Day 4-Pass Card and the Qube Mini 7-Day 4-Pass Multi-Med Card are filled and ready to use.